Control Plasma LMW Heparin

For In Vitro Diagnostic Use

Control Plasma LMW Heparin

Art. No. 82 34 92

CHROMOGENIX

INTENDED USE

For the quality control of chromogenic heparin assays.

REAGENT

A lyophilized powder of human plasma containing low molecular weight (LMW) heparin.

PACKAGE

Control Plasma LMW Heparin, Low	1 ml	4 vials
Control Plasma LMW Heparin, High	1 ml	4 vials

CAUTION: Each donor unit used in the preparation of human source reagent has been tested by FDA approved methods for the presence of Hepatitis B surface antigen and antibodies to HIV 1 and 2 and Hepatitis C and found to be negative. However, since no test can completely rule out the presence of these blood borne diseases, the handling and disposal of human source reagents from this product should be made with care. Handle as potentially infectious¹.

PROCEDURE

Reconstitute the content of each vial with 1.0 ml of NCCLS type II water² or deionized water, filtered through 0.22 µm. Allow to stand at room temperature for 20 minutes. Swirl gently before use.

STABILITY

Stability after reconstitution: 24 hours at 20-25°C, or 14 days at 2-8°C in the original vial.

EXPECTED VALUES

Each lot of the Control Plasma LMW Heparin is calibrated against the 1st International WHO Standard for LMW heparin using COAMATIC® Heparin.

Expected values - see separate sheet.

Each laboratory should establish its own mean and standard deviation and should establish quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to Westgard et al. for identification and resolution for out-of-control situations³.

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REFERENCES

1. RICHARDSON J H and BACKLEY W E. Eds. Biosafety in Microbiological and Biomedical Laboratories. US. Dept. of Health and Human Services, Public Health Service, HHS Publication No (CDC) 84-8395, Washington, D.C. 1984.
2. National Committee for Clinical Laboratory Standards. Specifications for reagent water used in the clinical laboratory, NCCLS Approved Standard: ASC-3.
3. WESTGARD J O, BARRY P L. Cost-effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.